

1. 產品規格

產品名稱	ECG 心率監測器 符合 ISO13485 標準
規格	3.5mm/4.0mm 雙孔 US 規格
長度	120cm 標準長度
重量	約 100g (含電池)
功能	TENS EMS 電療/ECG/EKG/EEG/EMG 監測
電池壽命	500h

US 規格產品-雙孔 ECG 心率監測器，符合 ISO13485 標準，適用於醫療、運動、康復等領域。

產品規格：雙孔規格，符合 US 標準。TENS, EMS, ECG 及 EMG 功能齊全，適用於多種電療及監測需求。

產品長度：標準長度為 120cm，可根據需求進行調整。3.5mm - 4.0mm 雙孔規格 — TENS 及 EMS 適用。

US 規格產品-雙孔 ECG 心率監測器，符合 ISO13485 標準，適用於醫療、運動、康復等領域。

產品規格：雙孔規格，符合 US 標準。TENS, EMS, ECG 及 EMG 功能齊全，適用於多種電療及監測需求。

產品規格：雙孔規格，符合 US 標準。TENS, EMS, ECG 及 EMG 功能齊全，適用於多種電療及監測需求。

產品規格：雙孔規格，符合 US 標準。TENS, EMS, ECG 及 EMG 功能齊全，適用於多種電療及監測需求。


產品規格：雙孔規格，符合 US 標準。TENS, EMS, ECG 及 EMG 功能齊全，適用於多種電療及監測需求。

2. 產品特點





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: □ □□□ □□ □□ □□□ □□□ □ □□□?

[illegible]

Q2: □□ □□ □□□ □□□ □ □□□□?

A: PayPal, T/T, ☐ ☐ ☐ L/C

Q3: 如何選擇合適的數據類型？

註: 本廠產品。通過 3C 認證。QC 認證。通過 CE, RoHS, SGS, FCC, ISO9001:2008 認證。

Q4: 何? 何! 何 何何 何何何 何何何 何 何 何何 何何 何 何何 何 何何?

A: 000000, 0000 1000 0000 0000 0000. 0000 OEM/ODM 0000 000000. 0000 0000 00 00 00 00 000000 0000 0000
0.

Q5: 如何 選擇 適合 的 數據 庫? 在 選擇 數據 庫 時 有 些 什麼 考慮 因素?

A: 2~4週間, 12週間 前後 程度. DHL, Fedex, UPS, TNT, EMS 利用 可能.

Q6: □□□ □□□ □□□□? □□□ □□□ □□□ □ □□□?

00: 00000. 0000 00 00 0000 00 0000 00000 00 00 0000. 00 0000 0 0000 00000 0000. 00 0000 0000 0000 0000 00
 100 80 A0