

Trade Term: FOB Shanghai/China

Quantity: 1,000,000 PC

Payment: T/T, L/C, Cash

Lead Time: 10-15 Days

Delivery: 3-7 Days, 15-20 Days

Certification: UL, CE, ROHS, FCC, UKCA


Shipping: Fedex, DHL, UPS

MOQ: 500PCS

Notes:

- [TENS 7000](#), Nicwell Care, TENS 3000 are available
- [TENS 7000](#): 10mm x 10mm x 10mm, 10mm x 10mm x 10mm.
- [TENS 7000](#) size: DC 2.35mm x 2.0mm x 1150mm. 1150mm.
- [TENS 7000](#) EMPI [Prospera](#) are available.
- [TENS 7000](#) are available!

Notes: [TENS 7000](#) are available:



00 00:

300 00 00 00 00 00 00 00

USB Type c - 2.0mm 0 00 00 00 00

Heymed Technology 000 0000 00000000. 000 000 0 0000, 0000 00 00 000 3.5mm 000 00 0000 00 00 000 00 000 00 000 00 00, 000, 00 00, 00 00 00 00 0000 00 000 00 000 0000.

00 000 1200 00 0 RoHS 0000 00 00 0 000 000000.

Hemed 00 ISO 13485 0000 00 00 00, ISO 9 0 0 1 00 00 00 00 0 ISO 1 4 0 0 1 00 00 000 00 0

□□□ UL, CQC, CE, SGS □□□ □□□□ □□□□□□.



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



📧📧📧 📧

Q1: 📧 📧📧📧 📧 📧 📧📧 📧📧 📧 📧📧?

A: 📧, Heymed📧 📧📧 📧 📧 📧📧 📧📧📧. 📧📧 📧📧📧 📧 📧📧 📧📧📧📧 📧📧. 📧📧 📧📧📧 📧 📧📧 📧📧 📧📧 📧📧📧.

Q2: 📧 📧 📧📧 📧📧 📧 📧📧📧?

A: PayPal, T/T, 📧 📧, 📧 📧 📧 L/C

Q3: 📧 📧📧 📧 📧📧 📧📧 📧📧📧 📧📧📧?

A: 📧📧 📧📧📧. 📧📧📧 📧 📧 📧📧📧 📧 QC 📧📧 📧📧📧 📧📧 📧 📧📧 📧 📧📧📧 📧📧📧. 📧 📧📧 Rohs, SGS, CE FCC, ISO9001:2008 📧📧 📧📧📧📧. 📧

Q4: 📧? 📧! 📧📧 📧📧📧📧 📧📧📧 📧 📧📧 📧📧 📧 📧📧 📧 📧📧?

A: 📧📧📧📧, 📧📧 100📧 📧📧 📧📧 📧 📧📧📧. 📧📧 OEM/ODM 📧📧📧 📧📧📧. 📧📧 📧📧 📧 📧 📧 📧 📧📧📧 📧📧 📧📧 📧.

Q5: 📧 📧📧 📧📧 📧📧📧? 📧 📧📧 📧 📧 📧📧📧 📧📧📧📧📧?

A: 📧📧 📧📧📧📧 📧 2-4📧📧📧 📧📧📧📧📧 12📧📧 📧📧 📧📧 📧 📧📧📧. 📧📧 DHL, Fedex, UPS, TNT, EMS📧 📧📧📧 📧📧📧📧.

Q6: 📧📧 📧📧 📧📧📧? 📧📧 📧📧 📧📧 📧 📧📧?

📧: 📧📧📧. 📧📧 📧 📧 📧📧📧 📧📧📧 📧 📧📧 📧📧📧 📧 📧 📧 📧📧 📧📧📧.