

Table Specifications

Connector	2mm Lemo 4P
Connector	3.5mm/4.0mm USB Type C
Length	120cm
Material	Medical grade plastic
Applications	TENS EMS ECG/EKG/EEG/EMG
Weight	500g

Notes:

LEMO B connector is not compatible with 2MM connector. Please use the correct connector for your device.

For more information, please visit [Lemo website](#), Ecg, Ekg, Eeg, Emg, Magnetic field, Magetic field.

Lemo connector is ECG compatible. [ECG EKG EMG EMC](#) compatible.

Heymed ISO 13485 ISO 9001 custom connector. Negative connector is not compatible.

Notes:





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

