

品名 規格


品名	00B4P lemo 2.0mm 4pin (0.5m)
規格	00B4P lemo 2.0mm 4pin
材質	PVC ,TPE,TPU
色別	黒, 白
長さ	0.5m, 1m
単位	個, 箱, 100個
MOQ	50PCS

備考:





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



歡迎來電

Q1: 貴公司主要經營哪些產品？

A: 我們, Heymed 主要經營各類醫療器械、美容儀器、家用電器、廚房電器、小電器等。我們擁有先進的生產設備和技術，產品質量優良，價格公道，歡迎各界人士垂詢。

Q2: 貴公司支持哪些付款方式？

A: PayPal, T/T, 信託, 銀行匯票 L/C

Q3: 貴公司產品通過哪些認證？

A: 我們通過了 CE, RoHS, SGS, CE FCC, ISO9001:2008 等認證。我們擁有完善的质量管理体系，確保產品質量符合國際標準。我們還提供 OEM/ODM 服務，歡迎來電諮詢。

Q4: 貴公司最小訂貨量是多少？

A: 我們支持 100 個起訂。我們提供 OEM/ODM 服務，歡迎來電諮詢。我們還提供代客組裝、包裝、貼標等服務，歡迎來電諮詢。

Q5: 貴公司產品交付時間是多久？

A: 我們產品交付時間為 2-4 周。我們擁有先進的生產設備和技術，產品質量優良，價格公道，歡迎各界人士垂詢。

Q6: 貴公司產品包裝規格是怎樣的？

A: 我們產品包裝規格為 100 個裝。我們提供 OEM/ODM 服務，歡迎來電諮詢。我們還提供代客組裝、包裝、貼標等服務，歡迎來電諮詢。