

Trade Term: FOB Shanghai/China

Quantity: 1,000,000 PC

Payment: T/T, L/C, Cash

Lead Time: 10-15 Days

Delivery: 3-7 Days, 15-20 Days

Certification: UL, CE, ROHS, FCC, UKCA

Shipping: By Air, By Sea, Fedex, DHL, UPS

MOQ: 500PCS

Notes:

- All items: New, Original, Used, Refurbished, etc. All items are guaranteed to be 100% functional and meet the specifications.
- Minimum Order: 500PCS. For bulk orders, please contact us for a special price.
- ED Items: New, Original, Used, Refurbished, etc. All items are guaranteed to be 100% functional and meet the specifications.
- Other Items: New, Original, Used, Refurbished, etc. All items are guaranteed to be 100% functional and meet the specifications. For more details, please contact us.

Heymed USB C - 4.0mm, lemo 2.0mm

Notes:



詳細情報:

Heymed Technologyは、医療機器の製造、開発、販売、およびサービスを提供する会社です。3.5mmのケーブルは、医療機器の接続に使用されます。このケーブルは、医療機器の接続に使用されます。このケーブルは、医療機器の接続に使用されます。

このケーブルは12ヶ月の保証期間があり、RoHSに準拠しています。

Hemedは、ISO 13485の認証を取得しており、ISO 9001の認証を取得しています。ISO 14001の認証も取得しています。また、UL、CQC、CE、SGSの認証も取得しています。



11111

Q1: □□ □□□□ □□ □□ □□□ □□□ □□□?

A: Yes, Heymed can be used for any kind of data. It can be used for any kind of data. It can be used for any kind of data.

Q2: □□ □□ □□□ □□□ □ □□□□?

A: PayPal, T/T, ☐ ☐, ☐ ☐ ☐ L/C

Q3: □□ □□□ □□□□ □□□ □□□ □□□ □□□□ □□□□?

A: 本公司產品。通過 ISO 9001:2008 認證。QC 部門嚴格控制產品質量。產品符合 RoHS, SGS, CE FCC, ISO9001:2008 標準。詳情請參閱說明書。

Q4: 如何? 为什么? 如何? 为什么? 如何? 为什么? 如何? 为什么? 如何? 为什么? 如何? 为什么?

A: 100% OEM/ODM available. We can provide you with the best price and quality.

Q5: 如何理解「物极必反，否极泰来」？

A: 2-4 weeks 12 weeks DHL, Fedex, UPS, TNT, EMS

Q6: 如何 選擇 合適 的 方法 來 解決 問題?

□: □□□□. □□□□ □□ □□ □□□□ □□□□ □□ □□□ □□□□□□.