





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

問と答

Q1: 貴社はどのような会社ですか？

A: はい。Heymedは、中国の主要な都市に拠点を設け、さまざまな種類の医療機器を製造する会社です。私たちは、高品質の医療機器を提供することを使命としています。

Q2: 貴社の主要な製品は何ですか？

A: はい、T/T、信用状、銀行振込などの支払い方法でL/C。

Q3: 貴社の製品はどのような規格に準拠していますか？

A: はい。貴社の製品は、QC、RoHS、SGS、CE、FCC、ISO9001:2008などの国際規格に準拠しています。

Q4: 貴社の MOQ (Minimum Order Quantity) は何ですか？

A: はい、100個以上です。OEM/ODMも受け付けています。お問い合わせください。

Q5: 貴社のリードタイムはどのくらいですか？

A: はい、2-4週間です。12週間以内です。DHL、Fedex、UPS、TNT、EMSなどの配送方法があります。

Q6: 貴社の製品はどのような保証がありますか？

A: はい。貴社の製品は、3年間の保証があります。1年間の8%のA1。