

EEG ECG

| | |
|-------------|--------------------------------------|
| Material | 100% Cotton, 100% Polyester ECG Lead |
| Color | Black, White, Blue Dupont Lead |
| Length | 30cm, 60cm, 90cm |
| Width | 1.5cm, 2.0cm, 2.5cm |
| Application | EEG, ECG, EMG, etc. |
| MOQ | 100PCS |

We are a professional manufacturer of ECG leads. We provide OEM ODM services. Our products are made of high-quality materials and are available in various colors and lengths. We also offer custom printing services. Our products are widely used in hospitals, clinics, and research laboratories. We are committed to providing high-quality products and excellent customer service. [DC 2.35](#) [EEG ECG](#) [USB](#).

Option A: Dupont Lead

Option B: 4.0mm Width Lead

Option C: [SCSI 20](#) Lead

Contact Us:





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

