

00 0000

000	00 00 00 00 000 000 00 EMS Adhesive 0 00
000	2.0mm 0 000
00	000
00	000
0000	00 00, 0000 00
00	00 00 000, 000 00 00
MOQ	100PCS

000 [2.0mm 0 0000 00 TENS 00 0000 00 00](#) 00 000 00 0000 00 0000 0000. 0000 00 00 0000 00 0 0000 00000 00 0 0000 00 00 000 0 00000. 0000 00000 TENS 00 00 00 00 3.5mm 0000 0000 00 00 00 00 0000 00000000.

000 0000 00 0000, 0000 0000 000000. Tens 00 00 00 0 00 0000 00 00 0000 00000 [00 Lemo 0000 0 0 00000 0 00 000000 00000 0](#), 00 00 0000, 00 0 00 000000.

00 00:





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

