

□□ □□: FOB □□/□□□/□□

□□ □□: □ 1,000,000 PC

□□: T/T, □□□□, □□ □□

Page 10 of 10

□□ □□: □□ □□ □□ 3-7□, □□ □□□ □ 15-20□

□□□ : UL, CE, ROHS, FCC, UKCA

□□ : □□, □□ □□ Fedex, DHL, UPS

MOQ:500PCS

11

1. 00 000 900, 00 000, 000 000 00 00 000 00 000 000 0 00 000 000 000 000000.
2. 00 00 000 00 000 0000 0000 0000 00 0000 000 0 0000. 0000 00 0 000 000000.
3. 000 00, 000 0000, TENS 00 0 00 00 00 000 000000. 0000 0000 0000 000 0 00 0000, 0000 000 0 0000 0000000.

111


90 DC2.35 1000 - 2.0 1000 ECG 1000, 1000 2.0 1000 1000 4 1000 1000

□□ □□ 12□□ □□ □ RoHS□ □□□□ □□ □□ □ □□□ □□□□.

Hemed® 符合 ISO 13485 质量管理体系, ISO 9001 质量管理体系及 ISO 14001 环境管理体系认证, 符合 UL, CQC, CE, SGS 等国际认证标准。



Figure 1 shows the dimensions of the microfluidic chip. The chip is a white, cylindrical device. The dimensions are indicated by dashed lines and arrows: the total length is 18.2mm, the diameter is 12mm, and the central channel diameter is 3.5mm.



A: 本公司產品。通過 ISO 9001:2008 認證。QC 部門嚴格控制產品質量。產品符合 RoHS, SGS, CE FCC, ISO9001:2008 認證。詳情請洽本公司。

