

	USB 3.0 Type C ~ 3.5mm
	USB C 3.5mm/4.0mm
	30cm
	TENS EMS / ECG/EKG/EEG/EMG
MOQ	100PCS

ISO13485

[ECG EKG EEG EMG](#) .

[USB C - ECG USB C 3.5mm ecg EEG ECG](#).
[\(TPU \) OEM ODM](#)),[ECG](#)
[M5](#) .


 , . [EEG ECG](#) . . .

 :





□□□□.□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



□□□□ □□

Q1: □ □□□□ □□ □□ □□□ □ □□□□?

A: □□□□. **Heymed** □ □□□ □□ □□ □□ □□□□, □□□ □□□ □□ □□□ □□□□ □□□. □□□ □□□□ □□ □□□ □□□ □□□□ □□□□ □.

Q2: □□ □□ □□□ □□□□□?

A: □□□, **T/T**, □□ □□, □□□ □□□ □ **L/C**

Q3: □□□ □□ □□□ □□□□ □□□ □□ □□□□?

A: □□□ □□□□□. □□□□ □□ □□ □□□□ □□ **QC** □□□ □□□□ □□□ □□ □□□ □□ □□□□ □□□□□. □□ □□□ **Rohs, SGS, CE FCC, ISO9001:2008** □□□ □□□□□□. □

Q4: □□? □! □□□ □□□ □□□□□ □□□□ □□□ □□□ □□□ □□□ □□□ □□□□?

A: □□□□□, □□□ **100** □ □□□ □□ □□□□□. □□□ **OEM/ODM** □□□□ □□□□□. □□□ □□□ □□ □□ □□ □□□□□ □□□□□□.

Q5: □□ □□ □ □□? □ □□□ □□ □□ □□□□ □□□□□?

A: □□ □□□ □ **2-4**□□□ □□ □ **12**□□ □□□ □□□ □ □□□□. □□□ **DHL, Fedex, UPS, TNT, EMS** □ □□□□ □□□□□.

Q6: □□□ □□□ □□□□? □□□ □□□ □□□ □ □□□□?

A: □□□□□. □□□ □□ □□ □□□□ □□ □□□ □□□□ □□ □□□□. □□ □□□ **3rd**□ □□□□. □□ □□□ □□□ □□□ □□□□ □□**1**□□ **8**□
A□ **1**□