

## 2.0 Electrode Needle Turn Dual 3.5mm 4.0mm Magnet Snap White Cable for Physical Therapy

illustrate:

Trade terms: FOB Shenzhen/Guangzhou/Hong Kong

Production capacity: 1,000,000 pieces per month

Payment method: T/T, Paypal, Western Union

Place of Origin: Dongguan, China

Delivery time: sample delivery time is 3-7 days, mass production is about 15-20 days

Certificate: UL, CE, ROHS, FCC, UKCA

Shipping: Air, Sea or Fedex, DHL, UPS

Minimum order quantity: 500PCS

Haimed Technology Magnetic Snap Button Factory provides high-quality magnetic snaps. It is not just a magnetic button. Haimed has always been a professional manufacturer of customized magnetic button wires in China, with supporting electrode lead buttons, for example, [3.5mm Stereo Jack 3-Lead ECG Snap-On Magnetic Electrode Cable](#), and [ECG Male Snap to 2.54mm DuPont Female Socket Connector Cable](#), and [20-pin connector to 4.0 mm electrode female ECG lead wire](#), and [USB A male to 3.5mm 4.0mm 5.0mm female snap-on electrode cable with switching power supply control](#) etc.

Product quality guaranteed for more than 12 months, compliant with RoHS and Reach standards and certified

Haimed has passed I S O 13485 medical device quality management certification, I S O 9 0 0 1 international quality management registration and I S O 1 4 0 0 1 environmental management system registration, and its products have passed UL, CQC, CE, SGS test reports