

Universal Adapter Wire 4.0mm Banana male plug to 10mm female Snap ECG Electrodes  
TENS unit Lead wires for electrical testing to illustrate:



Heymed Technology Test Leads Electrode Snap Button Factory Provide High Quality [Banana Plug to snaps cable adapter](#) Electrodes .

This customized [ECG /EKG electrode cable](#) with 10.0mm female snap is designed to avoid other external influences on the metal contact.

And it is widely used in medical industry,such as ECG EEG EKG EMG machine,we also do many similar items with DC3.5 plug to 4.0 electrode\_or with type c male and lemos connector cable and so on in custom.

What's more,Heymed support in custom of different colors,connectors,materials and so on.If any other requirements please feel free to contact us.



Product quality guaranteed for more than 12 months, compliant with RoHS and Reach standards and certified

Heymed has passed I S O 13485 medical device quality management certification, I S O 9 0 0 1 international quality management registration and I S O 1 4 0 0 1 environmental management system registration, and its products have passed the UL, CQC, CE and SGS test.