

Product name	2.5MM right angle audio male plug to two-in-one electrode ear clip from Guangdong professional customization ISO9001/13485, FDA, CE certification
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

#### application:

- **Hospital Physiotherapy Department:** Used in various medium and low frequency physical therapy devices to assist in the treatment of neuralgia, muscle strains and other diseases.
- **rehabilitation center:** Used in conjunction with rehabilitation equipment to stimulate patients' ear points and promote recovery.
- **home health care:** Can be connected to home physical therapy equipment, users can perform ear health massage at home to relieve fatigue.

#### Main features:

- **Unique Y-shaped splitter design:** The 2.5 audio elbow interface is firmly connected and not easy to loosen. The two-in-one ear clip design can perform physical therapy on both ears at the same time to improve efficiency.
  - **Customized services:** Guangdong's professional factory has strong strength and can customize products according to customer requirements, such as length, color, packaging, etc.
  - **Professional certification:** With ISO9001, ISO13485, FDA, CE (MDR 2017/745) and other certifications, the quality meets international standards and is safe and reliable.
- High-quality material: The wire is made of PVC material, which is soft, durable, anti-aging, corrosion-resistant, and adaptable to various usage environments.








### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

## **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer  
of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables,  
medical cables and chargers with good quality and best price.



## FAQ

### **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

### **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

### **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

### **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

### **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

### **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China