

Product name	2.5mm Audio Jack Dual Ear Clip Cable Reusable Electrode Lead Wire for Medical Health Devices China Dongguan Manufacture
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

Our 2.5mm audio jack dual ear clip cable is a high-quality electrode lead wire designed for health and medical electronic devices. It features a reliable 2.5mm male connector and reusable stainless steel ear clips with soft pads for patient comfort. Built with durable copper wiring and a flexible PVC jacket, this cable ensures long-lasting performance.

Manufactured by a professional factory based in Dongguan, China, we are an ISO13485 and ISO9001 certified supplier with rich experience in OEM and ODM production. This product is widely used in TENS therapy, EMS muscle stimulators, and other electrotherapy equipment. For similar applications, we also offer a [USB-C to 4.0mm female snap electrode cable with silicone earhook](#), which ensures reliable signal transmission and enhanced patient comfort, as well as a [custom 6.35mm audio plug to dual 2.0mm electrode pin cable](#) designed for professional physiotherapy settings.

Custom lengths, packaging, and private labeling are available upon request.





Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China