

Our 32-Channel EEG Cap is a high-performance, reusable solution compatible with the International 10-20 electrode placement system. Constructed from flexible and breathable materials, it ensures patient comfort during extended EEG recordings. The cap supports secure electrode placement and reliable signal transmission, making it suitable for medical, research, and educational use.. Made from high-quality elastic fabric, the cap offers a comfortable fit for long-term use and can be paired seamlessly with [high-quality EEG lead wires](#) to enhance signal transmission during monitoring sessions.

Manufactured by a leading EEG accessories supplier and factory based in Dongguan, China, we provide customization options to meet specific project needs, including different electrode materials, connector types, and sizes. To further improve the reliability of brainwave data collection, our cap can be used alongside [professional medical electrode patches](#), ensuring better contact and signal clarity. Reliable quality control and fast global delivery guarantee excellent user experience and performance. OEM and ODM services are available.

Key Features:

- Compatible with international 10-20 EEG standard
- 32 evenly spaced channels for high-resolution EEG signals
- Comfortable and stretchable fabric for long-term wear
- Reusable design with up to 50 wash cycles
- Customizable electrode and connector options
- Manufactured in Dongguan, China – OEM/ODM available






Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China