

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS


MOQ:500PCS

Features:

1. **High-Quality Construction:** Crafted from premium materials, this cable ensures durability and longevity, even in demanding medical environments.
2. **Precise Signal Transmission:** Engineered to deliver precise and reliable signal transmission, maintaining the integrity of medical data for accurate diagnostics and analysis.
3. **Universal Compatibility:** Compatible with a wide range of medical equipment, including EEG (Electroencephalography) machines, ECG (Electrocardiography) monitors, EMG (Electromyography) devices, TENS (Transcutaneous Electrical Nerve Stimulation) units, and more.
4. **Dual 2.0mm Pin Electrode Connections:** Features dual 2.0mm pin connectors, ensuring a secure and stable connection with medical electrodes for consistent and uninterrupted data acquisition.
5. **3.5mm Jack Interface:** Equipped with a standard 3.5mm jack interface, allowing seamless connectivity with devices such as computers, tablets, smartphones, and medical recording systems.
6. **Flexible Length Options:** Available in various lengths to accommodate different setups and configurations, providing flexibility in medical procedures and experiments.
7. **Easy to Use:** Plug-and-play design ensures hassle-free setup and operation, saving time and effort for medical professionals and researchers.

8. **Safe and Reliable:** Complies with industry standards for safety and performance, offering peace of mind during medical procedures and patient monitoring.
9. **Versatile Applications:** Suitable for a wide range of medical applications, including neurological studies, cardiac monitoring, muscle activity assessment, pain management, and rehabilitation therapy.
10. **Professional Grade:** Trusted by medical professionals, researchers, educators, and healthcare institutions worldwide for its quality and performance in clinical and research settings.

The [pin style electrode cable](#) with 3.5mm jack Drawing reference:



Next model: [2.0mm pin to 3.9mm male snap electrode adapter cable](#)

[USB Type c to 2.0mm pin style electrode cable](#)

Heymed Technology is a leading china manufacturer for [medical cables](#) and accessories, Magnetic snap [electrode cable](#) with 3.5mm audio jack, that is not only a certain specification of the magnet snap leads, Heymed always do the best for customization, which provide the different wiring specifications, cable, length custom, snap electrode leads number.

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Hemed has already certified by I S O 13485 Medical devices quality management, I S O 9 0 0 1 Registration for international quality management and I S O 1 4 0 0 1 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Office Team



Production technical team



R&D test team



Quality test team



Sample exhibition hall



Product warehouse

Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.