

3.5mm stereo jack 3/4 leads ecg snap magnetic electrode lead wire cable

Decription:

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China


Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

Drawing:



environmental management system and Products has Passed UL , CQC,CE, SGS Testing reports