


The 3.9mm Dry Electrode with Ag/AgCl Male Snap is a high-performance component designed for biomedical monitoring applications such as ECG, EEG, EMG, and neurostimulation. Constructed with a durable ABS housing and coated with Silver/Silver Chloride (Ag/AgCl), this electrode ensures excellent biocompatibility and stable signal acquisition without requiring conductive gels.

Its 3.9mm male snap interface provides secure and repeatable attachment to cables or electrode pads, making it suitable for both reusable and disposable applications—similar to our [Ag/AgCl Gel-Based ECG Electrodes](#), which offer reliable performance for short-term monitoring. Produced in our Dongguan-based factory, we offer customization options for color, packaging, and private labeling to meet OEM/ODM demands, just like our [Reusable Neurostimulation Electrodes](#) that support tailored designs for diverse therapeutic needs.

Features:

- Dry electrode – no gel required
- Durable ABS body with Ag/AgCl coating
- Male snap design for secure connectivity
- Consistent bio-signal transmission
- Made in China by a trusted supplier






Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China