

Product parameters

| | |
|------------------------|--|
| Product name | 4.0mm Medical Male Snap Silicone Sheath TO 4.0mm Medical Female Snap PVC Molding TENS EMS Devices Dongguan TENS EMS Lead Wire Manufacturer |
| connector | 4.0mm snap (male and female) |
| Cable length | 120cm or customized |
| Cable color | black or customized |
| Function | TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable |
| Minimum order quantity | 500pcs |

Manufactured by a professional ISO13485-certified TENS EMS lead wire supplier and factory in Dongguan, China, this custom cable features a 4.0 mm medical male snap with silicone sheath connected to a 4.0 mm female snap with durable PVC molding. Engineered for use with TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscle Stimulation) therapy devices, this cable ensures stable connection and high signal fidelity in therapeutic settings.

You might also consider our [2.0 mm Female Electrode Pin to 4.0 mm Female Snap Lead Wire for TENS EMS Devices](#), which provides a pin-to-snap format using medical-grade PVC or TPU, ideal for systems requiring connection to standard pin-type device outputs

The biocompatible silicone sheath offers flexibility and comfort during extended use, while the PVC-molded female snap end ensures a secure fit with most electrode pads. Ideal for OEM/ODM customization, this cable supports various wire lengths, snap spacing, and lead configurations.


Alternatively, check out our [Customized Plastic Lemo Plug to 4 Leads 4.0 mm Electrode Female Snap Cable](#), featuring a Lemo-compatible plug on one end and four 4.0 mm female snap connectors—designed for multi-lead applications like TENS, EMS, ECG, and EMG systems

Product photos:





Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China