

Product parameters

Product Name	6 Pin DIN Male Connector to Open Wire for Physical Therapy Machine Dongguan ISO13485-Certified Medical Cable Supplier
Connector	6 Pin DIN Male Connector
Cable Length	50cm or customized
Cable Color	White Black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
MOQ	500PCS


Manufactured by a professional ISO13485-certified medical cable factory in Dongguan, China, this 6 Pin DIN Male Connector to Open Wire Cable provides high-reliability performance for physical therapy machines, rehabilitation devices, and neuromuscular stimulation equipment.

As an example of our product range, please see our [Customized 6.35 mm Audio Male Plug to Dual 3.5 mm Female Snap Connector Cable with 2 Leads](#)—a high-quality therapy cable produced by our Dongguan ISO13485-certified plant for TENS/EMS applications.

The precision-engineered DIN male plug ensures secure electrical contact, while the open wire end allows for versatile connections and easy customization. Using medical-grade PVC or silicone insulation, the cable offers durability, biocompatibility, and superior flexibility, ideal for frequent clinical use.


As a trusted Dongguan medical cable supplier and plant, we provide custom solutions, including cable length, wire gauge, and connector configurations, meeting diverse OEM and ODM requirements for physiotherapy and rehabilitation applications. For example, we also manufacture the [4 Pin L-Shaped Metal LEMO Plug to 316 Stainless Steel Connector Therapy Cable](#), fully customizable and certified under ISO 13485, tailored for therapy devices such as TENS, EMS, and bioelectric stimulation systems.

Product Photo:





Company Info. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.