

3.5mm Audio Male Cable Silicone Earring to 4.0 Female Snap Tens Electrode Lead Cable For Data Transfer Task of The Physiotherapy Device

Product Picture :





Heymed Technology Test Leads Electrode Snap Button Factory Provide High Quality [Audio Male Cable straight Head to Electrodes Female Buckle](#).

Heymed always as a professional manufacturer of customized magnetic button wires in China, with supporting electrode lead buttons, such as [Right Angle 3.5mm Stereo Audio Jack To 3 Lead 2.5mm Ecg Snap Button ECG EEG Electrode Cable](#) ,Tens Electrode Lead Wire DC 2.35mm Male To Dual 3.5mm 4.0mm Ecg Snap Cable and [6.35MM TRS 1/4 Audio Jack To Dual 4.0MM Female ECG Snap Shiedling Wire Tens Cable](#) etc .

Heymed support in custom of different colors,connectors,materials and so on.If you have any quastions or need further details, please feel free to contact us.

Product quality guaranteed for more than 12 months, compliant with RoHS and Reach standards and certified

Heymed has passed I S O 13485 medical device quality management certification, I S O 9 0 0 1 international quality management registration and I S O 1 4 0 0 1 environmental management system registration, and its products have passed the UL, CQC, CE and SGS test.