

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

---

### Specification:

Ensure reliable and accurate electrode stimulator with our 2 leads 3.5mm 4.0mm [magnetic buckle lead wire](#) to 2 pole 6.35mm audio male wire cable. Designed for seamless connectivity and durability, this cable is an essential component for medical professionals and healthcare facilities.

### Features:

- **Compatibility:** Compatible with standard ECG equipment featuring 6.35mm mono jack audio male connectors.
- **Durability:** Constructed with high-quality materials for long-lasting performance in medical environments.
- **Ease of Use:** 3.5mm 4.0mm ecg magnetic buckle connectors securely attach to electrode leads for easy setup and removal.
- **Precision:** Provides clear and consistent power transmission for accurate TENS stimulation devices.
- **Length:** supported custom length suitable for various patient monitoring setups.

**Applications:** Ideal for hospitals, clinics, and medical practices requiring reliable ECG monitoring solutions. Suitable for use with a range of ECG machines and cardiac monitoring devices.

### Specifications:

- **Connector Type:** 3.5mm 4.0mm ECG magnetic buckle to 2 poles 6.35mm audio male
- **Number of Leads:** 2 leads
- **Cable Length:** 120cm Standard length (custom lengths available upon request)
- **Material:** High-quality, medical-grade materials

Heymed provide custom solutions with different connector for the electrode cable, Special [5-pole 3.5mm jack electrode lead wire](#), please contact with us.

---







Heymed Technology is a leading china ECG EKG cable manufacturer, Providing custom medical cable services, Heymed always do the best for customization, which provide the different wiring specifications, cable, length custom, snap electrode leads number.

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Hemed has already certified by I S O 13485 Medical devices quality management, I S O 9 0 0 1 Registration for international quality management and I S O 1 4 0 0 1 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing reports

---

## FAQ

**Q1: Can I start a sample order before I order more?**

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

**Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

**Q3: How is your factory doing regarding quality control?**

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

**Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?**

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

**Q5: Delivery time and warranty? What shipping service will you use for my products?**

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

**Q6: Where is your factory located? How can I visit there?**

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.