

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

Specification:

Description: Enhance your TENS (Transcutaneous Electrical Nerve Stimulation) and EMS (Electrical Muscle Stimulation) experience with this carbon conductive silicone rubber wristband. Designed for comfort and effectiveness, this wristband serves as an electrode, providing reliable electrical conductivity for therapeutic stimulation. Whether you're managing pain relief or muscle rehabilitation, this wristband ensures consistent and efficient electrical signal transmission.

Key Features:

- **Carbon Conductive Silicone:** Ensures optimal electrical conductivity and durability.
- **Wristband Design:** Comfortably wraps around the wrist for secure electrode placement.
- **Universal Compatibility:** Compatible with most TENS and EMS units using carbon electrodes.
- **Easy to Use:** Simply wrap around the wrist and connect to your device for immediate use.
- **Reusable and Washable:** Designed for multiple uses; easily cleaned with mild soap and water.

Specifications:

- **Material:** Carbon Conductive Silicone Rubber
- **Color:** Black (Other colors available upon request)
- **Size:** Adjustable to fit most wrist sizes

Applications:

- Ideal for pain relief therapy, muscle rehabilitation, and physical therapy sessions.
- Suitable for personal use or professional medical settings.

Notes:

Attached [Banana plug to 4.0mm ecg female snap electrode cable](#) sold seperately!



Heymed Technology is a leading china stimulator TENS EMS Electrode wristband manufacturer, Providing [custom medical cable](#) services, Heymed always do the best for customization, which provide the different wiring specifications, cable, length custom, snap electrode leads number.

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Hemed has already certified by I S O 13485 Medical devices quality management, I S O 9 0 0 1 Registration for international quality management and I S O 1 4 0 0 1 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing reports

FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.