

# Product parameters

Product Name	Cusomized 6 Pin Male Connector to 5 Leads 2.0mm Electrode Pin ECG Cable for Physical Therapy Machine   Dongguan Medical Cable Supplier
Connector	6 Pin DIN Male Connector
Cable Length	50cm or customized
Cable Color	White Black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
MOQ	500PCS

Manufactured by a professional ISO13485-certified medical cable supplier in Dongguan, China, this 6 Pin Male Connector to 5 Leads 2.0mm Electrode Pin ECG Cable is engineered for reliable performance in physical therapy, rehabilitation, and electrotherapy applications.

The precision-designed 6 Pin Male Connector ensures stable electrical contact with therapy machines, while the 5 lead wires terminate in standard 2.0mm electrode pins, providing quick and secure attachment to electrode pads. Using medical-grade PVC or silicone insulation, the cable delivers excellent flexibility, biocompatibility, and long service life. For related solutions, you may also explore our [DB15 Female Plug to 4.0mm Electrode Cable](#) designed for multi-lead therapy connections.

**Key benefits include:**

- **Stable Signal Transmission:** High-quality conductor materials reduce noise and signal loss.
- **Medical-Grade Safety:** Insulation materials meet international biocompatibility standards.
- **Versatile Applications:** Compatible with ECG monitoring, rehabilitation devices, and physical therapy equipment.
- **OEM & ODM Support:** Connector types, cable lengths, and colors can be customized according to client requirements.


As a trusted Dongguan medical cable manufacturer, we provide not only ECG cables but also a wide range of [tVNS Ear Clip Electrodes](#) and therapy accessories to meet different clinical and therapeutic needs.

## Product Photo:





## Company Info. Overview



## About Us

10+ Years Experienced Professional Manufacturer  
of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables,  
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

# FAQ

## **Q1: Can I start a sample order before I order more?**

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

## **Q3: How is your factory doing regarding quality control?**

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

## **Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?**

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use for my products?**

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

## **Q6: Where is your factory located? How can I visit there?**

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.