

Product parameters:

Product name	Customized 10 Pin pogo pin mangtic connector to 10 leads 4.0mm female snap ECG Cable
connector	3.5mm/4.0mm female and Pogo Pin Male
Cable length	80cm or customized
Cable color	white,grey,black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
Minimum order quantity	500pcs

Manufactured by a professional ISO 13485-certified medical cable supplier in Dongguan, China, this ECG cable combines precision engineering with robust performance. The pogo pin magnetic connector design ensures secure contact, low resistance, and consistent electrical transmission, reducing signal loss during monitoring.

The 4.0mm female snap leads are compatible with multiple electrode pads, ensuring versatility across a wide range of ECG and diagnostic equipment. Medical-grade PVC or silicone insulation provides flexibility and biocompatibility, ensuring safe and comfortable use during long-term patient monitoring. For multi-patient or multi-channel setups, see our [Multi Leads EEG ECG Monitoring Electrode Cable](#)

Looking for other compatible medical cables? Check our [9Pin Magnetic Spring Loaded Connector to 2.0mm Round Ring Terminal with 5 Lead Wire](#) for more options.

Features:


- Production name: 10Pin Pogo Pin to 4.0mm Electrode Female Snap 10 Leads Wire ECG Cable
- Connector 1: 10Pin Pogo Pin Male Connector
- Connector 2: 4.0mm ECG snap
- Cable color: black, white, gray
- Snap mold injection color: Red,Blue,Green, White.
- Heymed Located in Dongguan, [ISO 13485&ISO 9001 certified cusom cable assembly](#), we can produce ecg snap cable, nagetive plate and other cable assemblies, please click [Dongguan ECG Medical cable supplier](#) to learn more.

Product photos:





Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China