

Product parameters

Product Name	2 In 1 12P Din Connector Male to Test Clip Clamp Electrode Cable
Connector	12pin Din connector male,test clip clamp
Cable Length	30cm or customized
Cable Color	Medical gray or customized
Function	TENS EMS Transmission,ECG EEG ECG machine
MOQ	100PCS

This Circular DIN Connectors 12 Pin Male to Test Clip Clamp Electrode Cable is customized by our cilents,which is widely used in medical industry,with high quality meterial of medical grade shielded wire.Meanwhile,we have many kinds of [Custom Medical Cables](#),such as [Multi leads EEG ECG monitoring electrode cable](#),which with 8 channel active eeg electrodes cable and eeg passive cable and [ECG male snap to electrode grabber clip cable](#).


Any demands of this customized medical cable,please contact us.

Product Photo:





Company Info. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.