

Product parameters

Product Name	Customized 2 Pin Lemo Plug to 2.0mm Electrode pin Cable with 2 Lead Wires Dongguan ISO13485 Medical Cable Manufacturer
Connector	Plastic LEMO Connector,2.0mm Electrode pin
Cable Length	1m or customized
Cable Color	Black or customized
Applied equipments	EEG, ECG,EMG
Function	Transferring signal
MOQ	100PCS

Manufactured by an ISO 13485-certified medical cable supplier in Dongguan, China, this customized 2 Pin Lemo Plastic Plug to Dual 2.0mm Electrode Pin Cable is designed for precise and secure signal delivery in therapeutic and diagnostic devices like TENS, EMS, ECG, and EMG systems.

The Lemo-compatible plastic connector ensures a reliable fit with standard medical equipment, while the dual 2.0mm electrode pin terminals provide broad compatibility with electrode pads. Constructed with flexible TPU or durable PVC, the cable offers excellent tensile strength and patient comfort.

For users requiring audio-jack style connectivity, our [DC 2.35 mm Male Plug to Dual 3.5 mm Electrode Female Snap Cable](#) provides a compact and efficient solution for TENS/EMS devices, available in 30 cm standard length with custom options

Looking to integrate with modern USB-C devices? We also offer a [USB-C Male to 4.0 mm Female Snap Electrode Cable with Silicone Earhook](#), featuring a secure silicone earhook and durable construction — ideal for ECG/EMG monitoring setups

Whether you're selecting established options or planning fully customized assemblies, our Dongguan factory delivers consistent quality and performance backed by ISO 13485 certification.


.

Product Photo:





Company Info. Overview



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.