


Product Name	Customized 2.54mm 2 Pin Male Plug to 2.0mm Electrode Pin Lead Wire Medical ECG Cable
application	Stimulation machine, physical therapy tool
length	10 cm, 20 cm, 50 cm, 100 cm or custom
color	Black, gray, white, black and custom
Conductor material	Canned copper
pack	Packaged separately

Manufactured by a professional medical cable supplier in Dongguan, China, this 2.54 mm 2-Pin Male Plug to Dual 2.0 mm Electrode Pin Lead Cable is engineered for compatibility with various diagnostic and therapeutic devices such as ECG, TENS, EMS, and EMG systems. For applications requiring a slimmer footprint, our DC [2.35mm Male to DC 2.35mm Female Extension Cable](#) offers a more compact alternative while maintaining high signal fidelity. Likewise, those needing rugged durability can explore [4Pin 2.5mm pitch connector to 2.0mm electrode pin cable](#), built with reinforced strain-relief and thicker insulation for demanding clinical environments.

The 2.0 mm pin needle leads offer a secure interface with electrode pads, delivering consistent electrical signals during monitoring or therapy. Built with medical-grade materials such as TPU or PVC, the cable is both flexible and durable for repeated clinical use.

Custom lengths, colors, and connector configurations are available to meet specific OEM/ODM requirements. Ideal for hospitals, rehabilitation centers, and physiotherapy device manufacturers seeking reliable, high-quality electrode lead solutions.






Certification and compliance:

The cable meets stringent standards to ensure the highest level of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards of medical equipment quality management systems.
- **ISO 9001:** Make sure the product is manufactured under a sound quality management system.
- **FDA Approval:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Comply with the EU Medical Equipment Regulation (MDR) to ensure the safety and effectiveness of medical uses.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples, and if we have stock, custom samples will require installation costs. Thank you for asking us to check our quality for samples.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Question 3: How is your factory doing in quality control?

Answer: Quality is preferred. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has been certified by ROHS, SGS, CE FCC, ISO9001:2008. ETC

Question 4: Factory? Yes! Can I customize my favorite items, such as printing my own logo or modifying my design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Question 5: Delivery time and warranty? What transportation services will you use to ship my products?

A: Our delivery time is about 2-4 weeks and we can provide a 12-month delivery warranty. We use DHL, FedEx, UPS, TNT, EMS as customers.

Question 6: Where is your factory? How should I visit?

Answer: Of course. We are pleased to welcome you to Dongguan and show you our workplace. Our factory is located in the third one. A Building A, No. 8, located at Kangyang 1st Street, Shebei Village, Dongguacheng, Guangdong Province, China