

Product Name	Customized 2pin connector To gold-plated 2.0mm banana plug cable with 2 leads wire Dongguan TENS EMS Lead Wire Manufacturer
application	Stimulation machine, physical therapy tool
length	10 cm, 20 cm, 50 cm, 100 cm or custom
color	Black, gray, white, black and custom
Conductor material	Canned copper
pack	Packaged separately

Manufactured by a professional TENS EMS lead wire supplier in Dongguan, China, this 2-lead cable connects a 2-pin medical device output to dual 2.0mm gold-plated banana plugs, offering high conductivity and corrosion resistance for consistent electrotherapy performance.

The 2-pin connector ensures a stable and secure connection with most TENS and EMS units, while the gold-plated banana plugs provide excellent compatibility with standard electrode pads or stimulation accessories. The lead wire is jacketed with TPU or PVC medical-grade insulation, making it flexible, skin-friendly, and durable under frequent use.

This model shares similar build quality and application with our [DC 2.35mm male to dual 2.0mm electrode pin braided cable](#), which features a robust braided exterior for enhanced durability. Additionally, for clients needing alternative connector formats, our [2.54mm 2 pin to 2.0mm electrode lead wire](#) offers equally secure performance for TENS, EMS, and ECG systems.

Ideal for clinics, therapy centers, and device manufacturers, this cable can be customized in length, color, and packaging to meet OEM/ODM requirements.





Certification and compliance:

The cable meets stringent standards to ensure the highest level of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards of medical equipment quality management systems.
- **ISO 9001:** Make sure the product is manufactured under a sound quality management system.
- **FDA Approval:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Comply with the EU Medical Equipment Regulation (MDR) to ensure the safety and effectiveness of medical uses.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples, and if we have stock, custom samples will require installation costs. Thank you for asking us to check our quality for samples.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Question 3: How is your factory doing in quality control?

Answer: Quality is preferred. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has been certified by ROHS, SGS, CE FCC, ISO9001:2008. ETC

Question 4: Factory? Yes! Can I customize my favorite items, such as printing my own logo or modifying my design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Question 5: Delivery time and warranty? What transportation services will you use to ship my products?

A: Our delivery time is about 2-4 weeks and we can provide a 12-month delivery warranty. We use DHL, FedEx, UPS, TNT, EMS as customers.

Question 6: Where is your factory? How should I visit?

Answer: Of course. We are pleased to welcome you to Dongguan and show you our workplace. Our factory is located in the third one. A Building A, No. 8, located at Kangyang 1st Street, Shebei Village, Dongguacheng, Guangdong Province, China