

Product name	Customized 4.0mm Male Snap To Ear Clip Electrode Cable Dongguan Medical Cable Manufacturer
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black, Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

Manufactured by a professional medical cable supplier and factory in Dongguan, China. This 4.0 mm male snap to ear clip electrode cable is purpose-built for TENS, EMS, and microcurrent therapy systems. The cable features durable TPU or PVC insulation, offering both flexibility and strength in clinical or home environments.

The 4.0 mm snap connector securely interfaces with standard electrode pads, while the conductive silicone ear clips ensure gentle and effective signal transmission across the earlobes or auricular points—commonly used in auricular acupuncture or neurostimulation. For customers needing grabber-type variants, check out our [4.0 mm Electrode ECG Male Snap to Electrode Grabber Clip Cable](#), a CE/ROHS/FCC-certified, ISO13485-manufactured option with customizable cable lengths and shielding.

Additionally, we offer durable [Double-Sided Silicone Conductive Ear Clip Electrodes for CES Therapy](#), also ISO13485 certified, which work seamlessly with our snap-connector lead wires for effective auricular stimulation.

As an ISO13485-certified vendor, we support full customization in terms of cable length, connector type, packaging, and labeling, making this solution ideal for OEM/ODM partnerships in the medical and wellness device industries.

Application:

- **Hospital rehabilitation departments:** It is used in various rehabilitation treatment equipment to assist patients in muscle stimulation, nerve function restoration and other treatments.
- **Professional physiotherapy centers:** Provide convenient equipment connection solutions for physiotherapists, and support a variety of physiotherapy projects, such as electro - acupuncture therapy, electro - assisted tuina therapy, etc.
- **Home care:** Facilitate home users to carry out simple physiotherapy and health care on their own, such as relieving muscle pain and promoting blood circulation.
- **Research institutions:** As a connecting wire for experimental equipment, it is used for the collection and analysis of bio - electrical signals in medical research.





Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China