

Product name	Customized 90 Degree 3.5mm Audio Plug to 3 Leads 2.0mm Electrode Pin Cable for TENS/EMS Units ISO13485 Dongguan Manufacturer
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black,Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

The 3.5mm male audio plug features a right-angle design to reduce strain during use. It connects seamlessly to devices such as TENS units or EMS stimulators. On the other end, the three 2.0mm electrode pin leads are compatible with standard electrode pads used in therapy treatments.

Manufactured in Dongguan, China, this cable is produced in an ISO9001 and ISO13485 certified facility, ensuring medical-grade quality and consistency. Suitable for global OEM/ODM partnerships, we offer custom branding, cable lengths, and packaging solutions.

Whether you're a distributor, healthcare device supplier, or medical equipment factory, this cable provides reliable performance for electrotherapy systems. For applications requiring a dual-lead configuration, consider our [3.5mm Jack to Dual 2.0mm Pin Electrode Connection Lead Wire Cable](#), which offers similar high-quality construction and compatibility. Additionally, our [90 Degree Right Angle Audio Jack to 3.5mm Electrode Snap Cable for EMS Device](#) provides an alternative connection solution for EMS devices, featuring a right-angle design for ease of use.



.




Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China