

Product parameters

Product name	Customized DB15PIN Female Plug to SR Splitter to 4.0mm Female Snap Electrode Cable with 13 Leads wire Dongguan Medical Cable Supplier
connector	3.5mm/4.0mm female and DB 15PIN Female Plug
Cable length	120cm or customized
Cable color	white black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
Minimum order quantity	500pcs

Our DB15PIN Female Plug to SR Splitter to 4.0mm Female Snap Electrode Cable with 13 Leads provides a dependable solution for electrotherapy, rehabilitation, and diagnostic equipment. Designed with precision engineering, this cable ensures minimal signal loss and optimal electrical conductivity, critical for accurate monitoring and therapy.

Reliable Connection: The DB15 female plug and SR splitter ensure stable data and power transfer between devices and electrode pads. For similar multi-lead designs, you may also explore our [10 Pin Pogo Pin Magnetic Connector to 10 Leads 4.0mm Female Snap ECG Cable](#), engineered for high-performance diagnostic applications.

Flexible Application: Compatible with a wide range of ECG, EMG, TENS, and physiotherapy systems.

Durable & Safe: Medical-grade PVC or silicone insulation offers flexibility, biocompatibility, and long service life.

Customizable Design: Length, snap size, wire colors, and packaging can be tailored to customer needs. For additional electrode connection options, check our [Dupont 2.54mm Female Plug to 4.0mm Male Snap Electrode Cable](#), which combines versatility with secure electrical performance.


As a professional Dongguan medical cable manufacturer, we provide OEM/ODM solutions with strict quality control for hospitals, clinics, and medical equipment brands worldwide.

Product photos:





Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China