

Product name	Customized Dual Color ABS Silicone Ear Clip for 2.0mm Round Tube Medical Cable   China Dongguan Manufacturer
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

This ABS silicone ear clip is an ideal solution for securing 2.0mm round tube cables in audio, medical, and communication applications. With a dual-color design, it enhances both functionality and aesthetics. The combination of rigid ABS and soft silicone ensures durability, flexibility, and comfort for prolonged use.

Key Features:

- Premium material: ABS + silicone for strength and flexibility
- Perfect fit: Designed for 2.0mm round tube cables
- Customization available: Wire length, colors, and branding options
- Certified quality: ISO13485 & ISO9001 manufacturer
- Reliable supplier: Produced in Dongguan, China, with a strong supply chain

In addition to this ear clip, we also offer [medical-grade TPU cables with customizable connectors](#), ideal for various healthcare applications. For disposable solutions, our [EEG electrodes with soft silicone coating](#) provide enhanced comfort and performance in clinical settings.

For bulk orders and customization options, contact us today.




### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

### **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China