

# Product parameters

Product Name	Customized Dupont 2.54mm Female Plug to 4.0mm Medical Electrode Male Snap Cable   Dongguan ISO13485-Certified Medical Cable Supplier
Connector	3.5mm/4.0mm male snap to Dupont 2.54mm Female Plug
Cable Length	50cm or customized
Cable Color	White Black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
MOQ	500PCS

The Dupont 2.54 mm female plug ensures firm contact with diagnostic equipment, while the 4.0 mm medical electrode male snap provides a quick and stable interface with electrode pads. Using medical-grade PVC or silicone insulation, the cable remains flexible, biocompatible, and durable for extended use in clinical settings. For a different configuration with electrode accessories, refer to our [Gold-Plated 5.0 mm ECG Electrode Female Snap Buttons for Medical Connector](#) manufactured under the same quality standards.”


Manufactured by a professional ISO13485-certified medical cable plant in Dongguan, China, this cable offers high-performance electrical conductivity and secure connections for a variety of medical devices. For another variant of electrode lead wire, you can also check our [Customized 4.0mm Male Snap To Ear Clip Electrode Cable](#), also produced in our Dongguan facility

## Product Photo:





**Company Info. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



## FAQ

### **Q1: Can I start a sample order before I order more?**

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

### **Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

### **Q3: How is your factory doing regarding quality control?**

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

### **Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?**

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

### **Q5: Delivery time and warranty? What shipping service will you use for my products?**

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

### **Q6: Where is your factory located? How can I visit there?**

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.