

# Product parameters

Product Name	Customized Micro Male Plug To 3.5mm Magnetic Female Electrode Snap Cable   Durable ECG/EKG Medical Cable Supplier from Dongguan
Connector	3.5mm/4.0mm female snap to Micro 5P Male
Cable Length	50cm or customized
Cable Color	White Black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
MOQ	100PCS

Manufactured by an ISO 13485-certified medical cable supplier in Dongguan, China, this Customized Micro Male Plug to 3.5mm Magnetic Female Electrode Snap Cable is designed to deliver high-performance connectivity for ECG/EKG, TENS, and physiotherapy equipment. The micro male plug ensures a precise fit, while the magnetic 3.5mm female snap provides quick, secure attachment and stable electrical conductivity.

Constructed with medical-grade PVC or silicone insulation, this cable offers flexibility, biocompatibility, and resistance to frequent sterilization processes. The durable design extends service life, making it suitable for both hospital and home-use therapeutic devices. For users requiring alternative plug designs, we also provide the [Type-C Plug to SR Splitter to 10 Leads 4.0mm Female Snap Electrode Cable](#) for secure and reliable medical connections.


Additionally, our [Custom Silicone Sheath 4.0mm Male Snap to 4.0mm Injection Molded Female Snap Cable](#) offers another high-quality solution, designed to meet the demands of TENS and EMS devices with exceptional durability and conductivity.

## Product Photo:





## Company Info. Overview



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



# FAQ

## **Q1: Can I start a sample order before I order more?**

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

## **Q3: How is your factory doing regarding quality control?**

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

## **Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?**

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use for my products?**

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

## **Q6: Where is your factory located? How can I visit there?**

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.