

Product parameters

Product name	Customized Type-C Plug to SR Splitter to 10 Leads 4.0mm Female Snap Electrode Cable Dongguan Medical Cable Manufacturer
connector	3.5mm/4.0mm female and USB Type C Male
Cable length	120cm or customized
Cable color	white black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
Minimum order quantity	500pcs

Manufactured by a professional ISO13485 certified medical cable supplier and factory in Dongguan China , this Type C to SR to 10 Leads 4.0mm Female Snap Electrode Cable is engineered to meet the needs of modern medical and wellness technologies. In addition to our 10-lead configuration, we also offer our compact [USB Type-C to Dual 3.5 mm Magnetic Snap Electrode Medical Cable](#), a streamlined solution ideal for dual-channel stimulation setups , as well as the versatile [Type-C to 4.0 mm Female Snap Cable with 2 Leads](#), optimized for standard TENS/EMS electrodes and multi-site connections

The Type C male connector offers compatibility with advanced portable stimulation units, while the strain relief(SR) structure improves cable lifespan and wear resistance under repeated use. Each of the 10 leads is terminated with a 4.0mm female snap, allowing for precise, secure contact with standard electrode pads.


As a leading OEM/ODM medical cable manufacturer, we offer full customization on lead color, length, connector orientation, and material type (TPU, silicone, or PVC). The cable is especially suitable for multi-site therapy, biometric monitoring, or rehabilitation systems requiring multiple contact points.

Product photos:





Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China