

Product name	Customized USB Type-C to 4 Leads Crocodile Clip Cable Dongguan Crocodile Cable Manufacturer
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black, Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

Designed for precision signal transmission and flexible connectivity, this Type-C to 4 leads crocodile clips cable is ideal for testing, temporary power connections, wearable devices, and lab use. Each crocodile clip is color-coded and securely bonded to ensure firm grip on terminals, conductors, or electrodes.

To support your OEM/ODM needs, we offer our [USB C to 4.0 mm Female Snap Electrode Lead Cable with silicone ear-hook](#)—a premium example of our Type-C medical cable series. Manufactured in Dongguan, it features PVC + copper construction, 0.5 m (customizable) length, and ISO 13485-certified factory quality

We also recommend our [4-in-1 Audio Jack to 4-Lead ECG Snap Shielding Cable](#), showcasing our expertise in multi-lead signal cables. This 1 m TPU-insulated cable supports multiple electrode sizes (e.g. 3.5 mm, 4.0 mm snaps), is fully shielded, and produced in our Dongguan factory under strict ISO & ISO 13485 systems

Produced by an experienced Dongguan manufacturer, our cables feature high-quality materials such as pure copper conductors and reinforced insulation to ensure durability and signal integrity. The USB Type-C plug ensures broad compatibility with modern diagnostic, measurement, or medical devices.

Our Dongguan factory supports bulk and wholesale orders, providing customized lengths, colors, and clip types to meet varied application needs.

Application:

- **Hospital rehabilitation departments:** It is used in various rehabilitation treatment equipment to assist patients in muscle stimulation, nerve function restoration and other treatments.
- **Professional physiotherapy centers:** Provide convenient equipment connection solutions for physiotherapists, and support a variety of physiotherapy projects, such as electro - acupuncture therapy, electro - assisted tuina therapy, etc.
- **Home care:** Facilitate home users to carry out simple physiotherapy and health care on their own, such as relieving muscle pain and promoting blood circulation.
- **Research institutions:** As a connecting wire for experimental equipment, it is used for the collection and analysis of bio - electrical signals in medical research.






Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China