

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA


Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

Features:

- Heymed [Replacement Electrode lead Wires](#) Compatible with TENS 7000, Nicwell Care, TENS 3000
- High quality Material: Copper Foil premium wire Leads , Durable for long time using.
- This Electrode lead wire for TENS7000 : DC 2.35mm Safety plug, with 2.0mm Male pin type connector. cable length 1150mm.
- Our TENS lead wires DO NOT Fit EMPI or Prospera Devices
- Please feel free to contact us if you need any help !

The [China wholesale TENS unit replacement lead wire](#) TENS units replacement lead wire Drawing reference:



Heymed Technology is a leading china TENS electrode cable vendor for medical cables and accessories, Magnetic snap electrode cable with 3.5mm audio jack, that is not only a certain specification of the magnet snap leads, Heymed always do the best for customization, which provide the different wiring specifications, cable, length custom, snap electrode leads number.

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Hemed has already certified by I S O 13485 Medical devices quality management, I S O 9 0 0 1 Registration for international quality management and I S O 1 4 0 0 1 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Office Team



Production technical team



R&D test team



Quality test team



Sample exhibition hall



Product warehouse

Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.