


Product Name	DC2.35mm safety standard male plug to 2 leads 2.0mm electrode pin braided cable TENS EMS Lead Wire Manufacturer
application	Stimulation machine, physical therapy tool
length	10 cm, 20 cm, 50 cm, 100 cm or custom
color	Black, gray, white, black and custom
Conductor material	Canned copper
pack	Packaged separately

Manufactured by a professional medical lead wire supplier in Dongguan, China, this durable and safe DC 2.35mm to 2-lead 2.0mm electrode pin cable is built to meet international medical device standards. The safety-enhanced male DC plug ensures secure insertion into stimulation equipment, while the braided cable body provides extra flexibility and anti-tangle performance ideal for frequent clinical or home use.

Each [2.0mm pin connector](#) is compatible with widely-used electrode pads for electrotherapy systems. The cable is suitable for TENS, EMS, and similar devices, and is often used in pain management and physiotherapy clinics. Similar braided designs—such as our [2.0mm Electrode Needle to Magnetic Snap Cable](#)—are also available for varied connection interfaces.

We support full OEM/ODM customization, including private labeling, color coding, custom packaging, and cable length adjustments. All cables are RoHS-compliant and produced under ISO13485 certification, ensuring safe and stable signal delivery in demanding environments.






Certification and compliance:

The cable meets stringent standards to ensure the highest level of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards of medical equipment quality management systems.
- **ISO 9001:** Make sure the product is manufactured under a sound quality management system.
- **FDA Approval:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Comply with the EU Medical Equipment Regulation (MDR) to ensure the safety and effectiveness of medical uses.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples, and if we have stock, custom samples will require installation costs. Thank you for asking us to check our quality for samples.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Question 3: How is your factory doing in quality control?

Answer: Quality is preferred. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has been certified by ROHS, SGS, CE FCC, ISO9001:2008. ETC

Question 4: Factory? Yes! Can I customize my favorite items, such as printing my own logo or modifying my design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Question 5: Delivery time and warranty? What transportation services will you use to ship my products?

A: Our delivery time is about 2-4 weeks and we can provide a 12-month delivery warranty. We use DHL, FedEx, UPS, TNT, EMS as customers.

Question 6: Where is your factory? How should I visit?

Answer: Of course. We are pleased to welcome you to Dongguan and show you our workplace. Our factory is located in the third one. A Building A, No. 8, located at Kangyang 1st Street, Shebei Village, Dongguacheng, Guangdong Province, China