


Product Name	DIN 6 Pin to 3 Lead 2.0mm Electrode Pin Splitter Cable with Metal Protection   TENS EMS Lead Wire Manufacturer
application	Stimulation machine, physical therapy tool
length	10 cm, 20 cm, 50 cm, 100 cm or custom
color	Black, gray, white, black and custom
Conductor material	Canned copper
pack	Packaged separately

Manufactured by a professional medical cable supplier in Dongguan, China, this Din 6 Pin to 3 Lead 2.0mm Electrode Splitter Cable offers reliable performance and superior protection. The DIN plug is reinforced with a metal protective sleeve to withstand repeated insertion and harsh usage environments. For advanced connection requirements, you may also explore our [Lemo Connector to 2mm electrode lead wire with plastic protect cover](#), which supports complex signal transmission in professional medical systems.

The 3-way 2.0 mm pin output is compatible with standard electrode pads or lead wires used in TENS/EMS/EMG devices, ensuring accurate electrical stimulation. This cable is widely used in rehabilitation centers, physiotherapy clinics, and home therapy setups. For compact or wearable devices, consider our [6.35mm Audio Male Plug to Dual 3.5mm Female Snap Connector Cable](#), designed for lightweight and portable monitoring applications.





### **Certification and compliance:**

The cable meets stringent standards to ensure the highest level of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards of medical equipment quality management systems.
- **ISO 9001:** Make sure the product is manufactured under a sound quality management system.
- **FDA Approval:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Comply with the EU Medical Equipment Regulation (MDR) to ensure the safety and effectiveness of medical uses.

### **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples, and if we have stock, custom samples will require installation costs. Thank you for asking us to check our quality for samples.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

## **Question 3: How is your factory doing in quality control?**

Answer: Quality is preferred. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has been certified by ROHS, SGS, CE FCC, ISO9001:2008. ETC

## **Question 4: Factory? Yes! Can I customize my favorite items, such as printing my own logo or modifying my design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Question 5: Delivery time and warranty? What transportation services will you use to ship my products?**

A: Our delivery time is about 2-4 weeks and we can provide a 12-month delivery warranty. We use DHL, FedEx, UPS, TNT, EMS as customers.

## **Question 6: Where is your factory? How should I visit?**

Answer: Of course. We are pleased to welcome you to Dongguan and show you our workplace. Our factory is located in the third one. A Building A, No. 8, located at Kangyang 1st Street, Shebei Village, Dongguacheng, Guangdong Province, China