

Product name	Factory Custom 2.0 Tens Lead Wires to 2P Male Pin Physiotherapy Cable by ISO13485 Manufacturer Heymed - Manufacturer Direct
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black,Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

## Product Information:

Our [2.0mm Electrode Pin Extension Cable](#) provides reliable and stable connections for various testing and diagnostic applications. This high-quality cable features male-to-female 2.0mm electrode pins, which are compatible with a wide range of devices that require electrode connectivity. The nickel-plated brass electrode pins ensure excellent conductivity, while the [PVC insulated copper wire](#) provides flexibility and durability for long-term use.

## Key Features:

- **Male to Female Design:** Secure connection for extending electrode testing devices.
- **Durable Material:** PVC insulation and nickel-plated brass for enhanced conductivity and longevity.
- **Flexible and Reliable:** 1-meter cable length, with options for custom lengths to suit your needs.
- **Wide Compatibility:** Compatible with a variety of medical and laboratory electrode equipment.
- **High Current Capacity:** Supports up to 5A, ideal for medical and laboratory applications.

## Application:

- Medical diagnostics and electrotherapy
- Laboratory equipment and experiments
- Electrodes for ECG, EMG, EEG devices
- Electrical testing and measurement setups

## Why Choose Us?

As an experienced [Electrode Pin Cable Manufacturer](#) based in Dongguan, China, we specialize in producing high-quality, customizable cables for the medical, laboratory, and electronics industries. Our strict quality control process ensures every cable meets the highest standards of safety and performance.

For other solutions related to Electrode Testing Cables, Medical Device Connectors, and Electrode Extension Leads, check out our full range of products to meet all your testing and medical connection needs.






### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

## **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer  
of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables,  
medical cables and chargers with good quality and best price.



# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China