

| | |
|--------------------|--|
| Product name | Factory Custom 6.35mm Audio Plug to Dual 2.0mm Electrode Male Pin Physiotherapy Cable by ISO13485 Manufacturer Heymed - Manufacturer Direct |
| application | Stimulation machine, physical therapy instrument |
| length | 10cm, 20cm, 50cm, 100cm or customized |
| color | Black,Gray, white, black and customized |
| Conductor material | Tinned copper |
| Package | Individually packaged |

Application Scenarios□


- **Clinics:** Used by physiotherapists in professional clinical settings for patient treatment.
- **Home - based Rehabilitation:** Ideal for patients who are undergoing rehabilitation at home under the guidance of a medical professional.
- **Wellness Centers:** Employed in wellness centers for relaxation and minor pain - management therapies.

Main features:

- **Durable Connectors:** The [6.35mm audio plug cable](#) is made of high - quality metal, providing a secure and long - lasting connection to the audio source. It resists wear and tear from frequent plugging and unplugging.
- **Precise 2.0mm Electrode pin:** The two [2.0mm electrode pin cable](#) are accurately crafted to ensure efficient transmission of electrical impulses. Their size is optimized for precise placement on the body's treatment points.
- **Flexible Cable:** The cable in - between is flexible yet sturdy, allowing for easy movement during treatments without tangling or breaking. It is made of high - grade materials that can withstand bending and stretching.



.




Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China