

Product name	High Quality USB A to 4.0 Electrode Female Snap with Switch Physiotherapy Cable for TENS Machine Therapy Sessions
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black,Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

## Purpose and Application Scenarios

- **Physiotherapy Clinics** Ideal for professional physiotherapy clinics, it can connect various electrical stimulation devices with [Lemo 2P female plug cable](#) and JST 5P terminals, accurately delivering therapeutic currents.
- **Home Healthcare** Convenient for home - based physiotherapy. It connects the Lemo and [JST terminals cable](#) of home - use physiotherapy devices, ensuring continuous treatment and recovery for patients.
- **Rehabilitation Centers** In rehabilitation centers, it can be used to connect rehabilitation equipment with corresponding interfaces, assisting patients in regaining muscle strength, improving joint mobility, and relieving pain.

## Main features:

- Constructed with a Lemo plastic female connector that ensures a secure and reliable connection. The plastic material is durable, resistant to wear and tear, and can withstand frequent use in physiotherapy settings.
- Features a JST 5P terminal which is known for its precision and compatibility. The terminals are well - crafted to provide stable electrical conductivity, reducing the risk of signal loss or interruption.



.










### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

## **Company information. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10 years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China