

# Product parameters

Product name	Factory Custom USB C Male to Silver Chloride EEG Suction Cup and Ear Clip Lead Wires for TENS EMS Devices   Dongguan TENS EMS Lead Wire Manufacturer
connector	USB Type C Male
Cable length	120cm or customized
Cable color	white black or customized
Function	TENS EMS Transmission/Massager ECG/EKG/EEG/EMG Adapter Cable
Minimum order quantity	500pcs

Manufactured by an ISO13485-certified medical cable supplier and factory in Dongguan, China, this custom USB Type-C male to silver chloride EEG suction cup and ear clip cable is designed for high-fidelity bio-signal transmission in therapeutic and diagnostic devices.

The silver/silver chloride suction cups and ear clips ensure excellent electrode-skin conductivity and stable long-term contact. This configuration is ideal for EEG neuro-monitoring, electrotherapy, and TENS/EMS therapy applications. For similar configurations, our [Plastic Lemo Plug to 4 Leads 4.0mm Electrode Snap Cable](#) provides another option for TENS/EMS connectivity, while the 3.5mm Audio Plug to EEG Cup Electrode Cable is suitable for systems requiring standard audio jacks and silver chloride EEG electrodes.

The [USB-C connector](#) ensures broad compatibility with modern electro-medical devices, while the lead wire can be tailored in length, jacket material, and wire gauge to meet specific clinical or research needs.

## Product photos:










## Company information. Overview



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

# FAQ

## **Q1: Can I start a sample order before ordering more?**

A: Of course, Heymed provides free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, bank transfer, Western Union and L/C

## **Q3: How does your factory do in terms of quality control?**

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from the beginning to the end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## **Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?**

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use to ship my product?**

A: Our delivery time is about 2-4 weeks, and we can provide a 12-month warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

## **Q6: Where is your factory located? How can I visit there?**

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China