

Product parameters

Product Name	Factory Customized 4.0mm Banana Plug to 2.0mm Electrode Pin Cable High-Quality Connectivity Solution from Dongguan China
Connector	Banana 4.0mm plug,2.0mm Electorode Pin
Cable Length	1m or customized
Cable Color	Black or customized
Applied equipments	EEG, ECG,EMG
Function	Transferring signal
MOQ	100PCS

Our factory-customized 4.0mm banana plug to 2.0mm electrode pin cable is designed for reliable signal transmission in medical, therapy, and monitoring applications. Produced in Dongguan, China by an experienced cable manufacturer and supplier, this cable features high-conductivity copper wire and durable TPU or PVC insulation, ensuring both flexibility and longevity. The connectors are securely molded for strong, stable connection during use in TENS or EMS equipment. OEM and ODM orders are welcome, with full customization options for size, cable type, and packaging. We offer fast delivery and consistent quality as a direct medical cable vendor.

For applications requiring multi-lead configurations, our [Multi-Functional 18in1 USB Type C To 2.0 Electrode Pin ECG Lead Wire Cable](#) offers versatile connectivity options. Additionally, for extending the reach of TENS units, consider our [TENS Unit Replacement Lead Wire 2.0mm Pin Female Extension Electrode Cable](#), designed for seamless integration with existing systems.


As a Dongguan-based manufacturer, supplier, and factory, we are committed to delivering high-quality, customizable medical cables to meet diverse client needs.

Product Photo:





Company Info. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.