

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

Features:

- Whole Body Care: Treat multiple areas, including shoulders, elbows, wrists, lower back, thighs, knees, and even the feet. It stimulates growth factor release, repairs cells, and relieves joint and muscle pain.
- Safe and Effective: PSP10 delivers low-intensity shockwaves for a non-invasive, painless, anesthesia-free experience. Shockwaves trigger your body's self-healing, improve metabolism, and initiate natural recovery.
- ED Solution: Low-intensity shockwaves enhance local blood circulation, promote new blood vessels, and improve penile health, effectively treating erectile dysfunction.
- User-Friendly: With touch controls, it offers a one-touch medical-grade treatment experience. Quickly alleviate joint and muscle pain from prolonged sitting or exercise. Small in size, powerful in performance.

Heymed custom variety high quality [medical cables](#) in China, we provide [USB C to 4.0mm snap electrode cable](#), [lemo connector to 2.0mm pin electrode cable](#) Etc.

Product Photo:



Company Information:

Heymed Technology is a leading China manufacturer for medical cables and accessories, Magnetic snap electrode cable with 3.5mm audio jack, that is not only a certain specification of the magnet snap leads, Heymed always do the best for customization, which provide the different wiring specifications, cable, length custom, snap electrode leads number.

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Hemed has already certified by I S O 13485 Medical devices quality management, I S O 9 0

0 1 Registration for international quality management and I S O 1 4 0 0 1 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing reports.



FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC,ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.