

| | |
|--------------------|--|
| Product name | Guangdong Factory Custom 3.5MM Audio Male Connector to Alligator Clips Physiotherapy Wire for Therapy |
| application | Stimulation machine, physical therapy instrument |
| length | 10cm, 20cm, 50cm, 100cm or customized |
| color | Black, Gray, white, black and customized |
| Conductor material | Tinned copper |
| Package | Individually packaged |

Main features:

- **Interface Advantages:**The 3.5MM audio male connector has strong versatility and is compatible with a variety of signal sources. The [alligator clips therapy cable](#) have a firm grip and are easy to operate.
- **Flexible Customization:**Cable length, color, and packaging can be customized according to customer requirements to meet diverse procurement needs.
- **Safe and Durable:**Environmentally friendly insulating materials are selected to prevent electric leakage, and high - quality copper cores ensure stable signal transmission. The structural design of the [3.5mm audio male](#) connector has been optimized, and the internal contact points are made of highly conductive materials. When combined with high - quality copper - core cables, it can minimize signal transmission loss.

Application:

- **Family Use:** It can be easily connected to household massagers and small - sized acupuncture and moxibustion devices, meeting the daily health - care physiotherapy needs.
- **Small - sized Clinics:** It is compatible with various instruments, and its customizable features help clinics provide personalized treatments.
- **Rehabilitation Centers:** With its stable performance, it can adapt to the high - frequency use in rehabilitation centers and ensure the therapeutic effects.










Certification and Compliance:

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- **ISO 9001:** Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China