

Product Name	M12 M16 4 Pin Male Connector to Dual 4.0mm Magnetic Female Snap Lead Wire Cable for ECG ISO13485 manufacturer Dongguan
application	Stimulation machine, physical therapy tool
length	10 cm, 20 cm, 50 cm, 100 cm or custom
color	Black, gray, white, black and custom
Conductor material	Canned copper
pack	Packaged separately

Manufactured by a professional medical cable supplier in Dongguan, China, this 4-Pin M12 or M16 connector cable features dual 4.0 mm magnetic female snap terminals, making it ideal for physiotherapy and rehabilitation systems. Similar magnetic snap cables, such as our [2.0 Electrode Needle to Dual 3.5 mm/4.0 mm Magnetic Snap Cable](#), are also widely used in TENS and EMS applications due to their compatibility and ease of use.

The strong magnetic design ensures quick connection with electrode pads while preventing accidental detachment. The cable's TPU or PVC jacket is resistant to flex fatigue and ensures safe use in clinical or home environments. The M12/M16 4-pin plug provides a secure interface with many medical stimulation devices, including those used for TENS, EMS, ECG, and EMG applications, just like our [Lemo Plug to 4.0 mm Gold-Plated Snap Lead Cable](#), which is commonly selected for ECG and rehabilitation systems.

Fully customizable in length, color, and packaging, this product is manufactured in compliance with ISO 13485 and RoHS standards, supporting OEM and ODM services.





Certification and compliance:

The cable meets stringent standards to ensure the highest level of quality and safety in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards of medical equipment quality management systems.
- **ISO 9001:** Make sure the product is manufactured under a sound quality management system.
- **FDA Approval:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Comply with the EU Medical Equipment Regulation (MDR) to ensure the safety and effectiveness of medical uses.

Company information. Overview



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

FAQ

Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples, and if we have stock, custom samples will require installation costs. Thank you for asking us to check our quality for samples.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Question 3: How is your factory doing in quality control?

Answer: Quality is preferred. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has been certified by ROHS, SGS, CE FCC, ISO9001:2008. ETC

Question 4: Factory? Yes! Can I customize my favorite items, such as printing my own logo or modifying my design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

Question 5: Delivery time and warranty? What transportation services will you use to ship my products?

A: Our delivery time is about 2-4 weeks and we can provide a 12-month delivery warranty. We use DHL, FedEx, UPS, TNT, EMS as customers.

Question 6: Where is your factory? How should I visit?

Answer: Of course. We are pleased to welcome you to Dongguan and show you our workplace. Our factory is located in the third one. A Building A, No. 8, located at Kangyang 1st Street, Shebei Village, Dongguacheng, Guangdong Province, China