

# Product parameters

Product Name	M16 M12 4Pin Connectors Plug Cord Wire Nurse Call Button Waterproof Cable for Patient   ISO13485 & ISO9001 Certified   Dongguan Manufacturer & Supplier
Connector	Plastic LEMO Connector,4.0mm ECG Female Snap Button
Cable Length	1m or customized
Cable Color	Black or customized
Applied equipments	EEG, ECG,EMG
Function	Transferring signal
MOQ	100PCS

Our M16 M12 4Pin Connectors Plug Cord Wire is designed to provide seamless connectivity for nurse call button systems, enabling patients to quickly and reliably signal for assistance. The 4Pin M12 connectors ensure secure connections, minimizing the risk of signal interruptions, while the M16 connectors offer easy installation and compatibility with a wide range of medical devices. For comparable solutions in medical connectivity, explore our [Integrated multi-lead electrode cable](#), engineered for precision in multi-channel monitoring systems with customizable configurations.

With an IP67 waterproof rating, this cable is resistant to water immersion and can withstand regular cleaning and disinfection procedures, making it suitable for use in healthcare settings where hygiene is critical. The medical-grade PVC jacket is non-toxic, phthalate-free, and easy to clean, meeting the strictest healthcare safety standards. For environments requiring durable and hygienic wiring, consider our [Custom Medical EKG Cables](#), designed to deliver reliable performance in diagnostic equipment with ISO 13485-certified quality assurance.

## Product Photo:










## **Company Info. Overview**



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years  
of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved  
manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



# FAQ

## **Q1: Can I start a sample order before I order more?**

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

## **Q2: What payment methods do you accept?**

A: PayPal, T/T, Bank Transfer, Western Union and L/C

## **Q3: How is your factory doing regarding quality control?**

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

## **Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?**

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

## **Q5: Delivery time and warranty? What shipping service will you use for my products?**

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

## **Q6: Where is your factory located? How can I visit there?**

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace. Our factory is located in 3rd. Floor, Building A, No.8 Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong, China.