

Trade Terms: FOB shenzhen/Guangzhou/Hongkong

Production Capacity:1,000,000 pcs per month

Payment : T/T, Paypal, western union

Place of Origin: Dongguan, China

Delivery time: sample lead time 3-7days, mass production is around 15-20 days

Certificates□UL, CE, ROHS, FCC, UKCA

Shipping: By air, by sea or Fedex, DHL, UPS

MOQ:500PCS

This is Negative Electrode Pads of TENS electrodes can be used multiple times for optimal relief from discomfort. Simply use the resealable package to keep your pads fresh and clean.

We have a variety of models and sizes available to suit the needs of different patients and surgical sites.

And We can supply OEM service to you ,make your LOGO's label to build your brand .

We also have some others electrode pad for physiotherapy like :[Silicone Electrosurgical Negative Plate With 6.3mm Plug](#) And [Lemo Plug Electrodes Pads Plate](#).

Products Quality Guaranteed 12months+ and RoHS and Reach Compliant & Certification

Heymed has already certified by ISO13485 Medical devices quality management, ISO9001 Registration for international quality management and ISO14001 Registration for environmental management system and Products has Passed UL , CQC,CE, SGS Testing reports.





FAQ

Q1: Can I start a sample order before I order more?

A: Sure, Heymed provide free sample if we have stock, customized sample need to be charged setup cost. We appreciate you asking us for a sample to check our quality.

Q2: What payment methods do you accept?

A: PayPal, T/T, Bank Transfer, Western Union and L/C

Q3: How is your factory doing regarding quality control?

A: Quality is priority. QC department with professional quality assurance specialist attaches great importance to quality control from the very beginning to the very end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 authentication. Etc

Q4: Factory? Yes! Can I customize the items I like, such as print our own logo or modify the design?

A: Yes, we are a factory with more than 100 employees. We offer OEM/ODM service. Please contact us or our sales representative for details.

Q5: Delivery time and warranty? What shipping service will you use for my products?

A: Our lead time is around 2-4 weeks and we can offer 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as the customer.

Q6: Where is your factory located? How can I visit there?

A: Sure. We are more than happy to welcome you to Dongguan and show you around our workplace.